

**Continuous dual-isocenter imaging:  
simultaneously extending the longitudinal  
and lateral intraoperative 3D CBCT field-of-  
view for assessing musculoskeletal trauma**

T. Reynolds<sup>1</sup>, S. Hatamikia<sup>2,3</sup>, O. Dillon<sup>1</sup>, Y. Ma<sup>4</sup>, A. Kanawati<sup>5</sup>,  
J.W. Stayman<sup>4</sup> and R. O'Brien<sup>1</sup>

2022 AAPM Annual Meeting – Washington, D.C.

Email: [tess.reynolds@sydney.edu.au](mailto:tess.reynolds@sydney.edu.au)

Twitter: @drtesreynolds

<sup>1</sup>University of Sydney, Sydney, Australia.  
<sup>2</sup>Austrian Center for Medical Innovation and Technology (ACMIT), Wiener Neustadt, Austria.  
<sup>3</sup>Medical Imaging Analysis and Artificial Intelligence (MIAAI) group, Donau-Private University Krems, Austria.  
<sup>4</sup>Johns Hopkins University, Baltimore, USA.  
<sup>5</sup>Westmead Hospital, Sydney, Australia.



**Continuous Dual Isocenter**



**Key Points and Conclusions**

Single volume  
longitudinal and lateral coverage

**95 × 44 cm**

Compared to stitched conventional  
acquisitions:

**80%** Reduction in scan time

**50%** Reduction in total  
projections acquired

ACRF IMAGE INSTITUTE

Email: [tess.reynolds@sydney.edu.au](mailto:tess.reynolds@sydney.edu.au)

Twitter: @drtesreynolds

